

sempermed® 
classic

CLASSIC THE ESTABLISHED CLASSIC

IDEAL TACTILE PROPERTIES –
specially developed dipping process
ENHANCED SAFETY – ideal compounding
EXCELLENT GRIP – high comfort

100+ 
YEARS
- EST. 1920 -
sempermed®
YOUR GLOBAL GLOVE EXPERT™

SEMPERMED® CLASSIC

THE ESTABLISHED CLASSIC

IDEAL TACTILE PROPERTIES

The specially developed dipping process guarantees the same wall thickness from fingers to cuff. This results in ideal tactile sensitivity and less hand fatigue even during longer operations.

SKIN CARE

The glove is powdered with pure cornstarch which makes it easy to don. The multiple washing process (leaching) makes the natural latex glove skin friendly and very comfortable to wear.

INCREASED SAFETY

Sempermed surgical gloves are classified as Medical Device (MD) as well as Personal Protective Equipment (PPE). This means that the glove can also be used for handling chemicals*.

EXCELLENT GRIP AND COMFORT

Sensitivity and a safe grip are essential for a successful surgery. This glove has a micro-rough non-slip surface to provide a good grip when handling instruments.



TECHNICAL DATA

sempermed® classic

TYPE	Sterile surgical glove for single use, powdered		
MATERIAL	Natural rubber latex (NR)		
COLOUR	Natural white		
GLOVE SHAPE	Fully anatomical with rolled rim		
SIZE/OVERALL LENGTH as per EN 455-2	5 ½, 6 and 6 ½ 7, 7 ½ and 8 8 ½ and 9	median 260 mm median 270 mm median 280 mm	
WALL THICKNESS measured in single layer in the palm area	min. 0.14 mm		
BARRIER PERFORMANCE as per EN 455-1	AQL 0.65		
FORCE AT BREAK as per EN 455-2	≥ 9 N		
DURABILITY in original packaging if stored as per DIN 7716, ISO 2230	5 years		
STERILISATION	Radiation STERILE R		
PACKAGING	Left and right glove with turned up cuff in inner pouch, sealed in medical peel pack. In dispenser carton with sterilisation indicator: 70 pairs In transport carton with sterilisation indicator: 420 pairs		
LABELLING	EN 1041, EN 455, EN 420, EN 374, ISO 15223		
PURPOSE according to MP Regulation 93/42/EWG according to PPE Regulation (EU) 2016/425	Medical device class IIa Single-use protective glove – PPE category III (protection against chemical substances for limited time)*		
ARTICLE NUMBERS	Size 5 ½ 826054520 Size 6 826054600 Size 6 ½ 826054620 Size 7 826054700	Size 7 ½ 826054720 Size 8 826054800 Size 8 ½ 826054820 Size 9 826054900	



CONTACT US!

Semperit Investments Asia Pte. Ltd.
8 Jurong Town Hall Road, #29-03/06 The JTC Summit,
Singapore 609434, Singapore
Tel.: +65 6274 4861 · Fax.: +65 6274 6977

Sempermed Europe GmbH
Am Belvedere 10,
A-1100 Vienna, Austria
Tel.: +43 1 79 777-0 · Fax: +43 1 79 777-630

sempermed@semperitgroup.com · www.sempermed.com

* Please note the Sempermed Chemical Resistance list is available on www.sempermed.com/chemical-resistance-brands

IMPORTANT NOTE: The latest product information is available at www.sempermed.com under the «Products» tab. Failure to observe this information, in particular with regard to (chemical) resistance, frequency of use and tolerability of the gloves, can result in personal injury and/or material damage. **Semperit is not liable for incorrect use of the gloves nor will Semperit be liable for improper storage and handling. In case of damage or unusual signs of wear and tear, gloves shall not be used.** In case of any doubt, expert advice should be sought before use. The information and classification is up to date as of the date of publication. Technical details are average values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors and amendments. **CAUTION: Natural latex can cause allergic reactions, including anaphylactic shock. Semperit bears no liability for allergic reactions.** ©Copyright 2022 Sempermed